

MIS
BIOCOM
SEVEN

Mount-less Implants



Overview.

MIS is proud to announce, that all MIS BIOCOM and SEVEN implants will now be supplied without mounts. These successful implants remain unchanged, but will be supplied without the blue mounts, and will be packaged in distinctive and easily recognizable boxes. These boxes will include a pull tab for quick and easy opening, and a new titanium sleeve, inside the tube, for maximum sterilization. This change will improve the handling, speed and safety of our implants. The BIOCOM and SEVEN mount-less implants, have their own tools, available as a set (MK-0037) or purchased separately.



Delivery tools for BIOCOM mount-less and SEVEN mount-less implants.

Delivery Tools.



MT-MRH20
Hex. Ratchet Long Adapter
for Int. Hex. Connection




MT-MRH10
Hex. Ratchet Short Adapter
for Int. Hex. Connection



MT-MMH20
Long Motor Adapter for
Int. Hex. Connection




MT-MMH10
Short Motor Adapter
for Int. Hex. Connection



MT-HHK24
Hex. Hand Key for
Int. Hex. Connection



MT-DI001
Implant Direction
Indicator



Angulation Indicator Probe



The logo for MIS (MIS Implants Technologies) features the lowercase letters 'mis' in a stylized, modern font. The 'i' has a vertical line extending upwards from its top, and the 's' has a long horizontal tail that curves slightly.

MIS Implants Technologies Ltd
www.mis-implants.com



0483 | ISO 9001:2000
ISO 13485:2003

Authorized European Representative:

MIS Implants Technologies GmbH
Paulinenstr. 12a, 32427 Minden, GERMANY
Tel: +49 571 97 27 69 0; Email: service@mis-implants.de