

Orthogum

Special orthodontic, very high viscosity C-Silicone.
Putty consistency.



Characteristics

- Fast setting
- Very high elasticity
- High tear strength
- Yellow colour
- Pleasant flavour with emetic reflex reduction

Technical Data

ADA 19	Type 1 Very High Viscosity
ISO 4823	Type 0 Putty Consistency
Mixing time	30"
Total working time (23 C/73° F)	1'
Time in mouth	1'
Setting time (23° C/73° F)	2' 30"
Compressive deformation(min-max)	2-5%
Elastic recovery	>99%
Dimensional stability	<-0.35%

Advantages

- Precision impressions with fixed orthodontic appliances
- Extremely easy and fast mixing with gel catalyst
- Good dimensional stability and possibility to develop the model even after 48 hours
- Possibility of obtaining multiple models with same Impression
- Safe disinfection with Zhermack Sterigum

Code

C100550 - 1 tub, 1500ml (2.46kg) MRP. 1800/-

C100551 - 1 tub, 1500ml + 2 Indurent gel .. MRP. 2300/-