




Pulp Master

- Check the vitality the pulse stimulus that is supplied to tooth.
- Press and hold speed increment key
- Display of pulse reading corresponding to maximum stimulus to 50 Indicates tooth condition
- Stimulus below 20 tooth is inflamed
- Stimulus between 21 to 35 indicates healthy tooth
- no stimulus seen at 35 indicates dead tooth



RELIABLE TECHNOLOGY

CONFIDENT DENTAL EQUIPMENTS LTD.

Pete channappa Indl.Este, Kamakshipalya, Magadi Road, Bangalore - 560 079

Ph : +91 80 23280801, 23280828 FAX : +91 80 23280822, Toll Free : 1800-425-3366.

marketingmgr@confidental.org, www.confidental.org