



Matrix retainer

Tofflemire retainer: Used to retain the matrix band in position during restoration of the Class II cavities or in case of Class I cavities with buccal or palatal extension.



TDV - Individual metallic matrix (with its own fixing system):

They are pre-contoured and have an incorporated fixing system that allows the regulation of the matrix according to the tooth size, fixing it firmly and dispensing the use of a matrix holder.



TDV - Individual plastic matrix (with its own fixing system):

They are pre-contoured and have an incorporated fixing system that allows the regulation of the matrix according to the tooth size, fixing it firmly and dispensing the use of a matrix holder.





CONFIDENT
SALES INDIA PVT LTD

Unimatrix:

- Due their convexity and slim thinness of 0,05mm they allow to rebuild perfect and tight contact points
- Three different sizes for a better adaptation to all clinical situations
- The rings can be sterilized by any method



Cervical matrix

Specially designed to condense and shape composites in class V restorations. Its form and flexibility can be easily adapted to the anatomy of the cervical area of the tooth when compressed against it.





CONFIDENT
SALES INDIA PVT LTD

Matrix Bands:

TDV Matrix retainer no 1 :



TDV Matrix retainer for molars :



TDV Matrix retainer for premolars:



Their high ductility and a thickness of 0,05mm allow these matrices to be contoured and adapted to the tooth very easily

Matrix strips/rolls;

TDV Metallic matrix strips(rolls):

- Dead soft stainless steel strips, with a thickness of 0,05mm
- Available in rolls of 3mts and 5mts length, in widths of 5mm, 6mm and 7mm



TDV Polyester matrix strips(rolls):

- Indicated for protecting adjacent teeth while the acid conditioning is being done.
- Due its transparency, resins can be cured through it



CONSERVATIVE ENDODONTICS



CONFIDENT
SALES INDIA PVT LTD

Wedges:

Anatomic wedges

- Anatomically shaped for a better adjustment to the interproximal spaces and to separate adjacent teeth
- Available in 5 sizes codified by colors



Reflexive wedges (transparent wedges)

- Cure through wedges that in combination with transparent matrices allows light curing composites on interdental cervical areas
- Available in 2 sizes



Elastic wedges

- Compatible with all matrix systems for Class II restorations
- Act simultaneously at vestibular and lingual regions, providing a perfect adjustment of the matrix to the tooth
- Promote a teeth separation without harming the gum, allowing to get a tight interproximal contact point
- Available in three sizes codified by colors



CONSERVATIVE & ENDODONTICS



CONFIDENT
SALES INDIA PVT LTD

Finishing & polishing strips:

Abrasive steel strips:

- For amalgam finishing
- Only 0,10mm thickness
- Available in widths of 4mm and 6mm
- Autoclavable

Diamond strips

- Indicated for orthodontic stripping / slices, and the finishing of filling materials
- Available in 2,5mm width
- Autoclavable

Finishing sandpaper strips

- For finishing and polishing composites
- Medium and fine grits separated by a neutral space for an easier insertion.
- Available in 2,5mm and 4,0 mm widths





Polishing paste:

Poligloss

Poligloss provides an excellent gloss on resins, amalgam, and metals. It is an aluminum oxide based polishing paste for general use in odontologic clinics and laboratory works

- Extra-fine grit (1 to 2 microns)
- Hydro-soluble & Easy cleaning



Diamond gloss

Diamond polishing paste was designed for produce a high gloss lustrous finish on all types of restorative materials as composites, ceramics, gold, etc.

- Useful at the clinic and the laboratory
- Extra-fine grit (1 to 2 microns)
- Hydrosoluble, easy cleaning





Finishing & polishing agents:

Optimize

- One single tip is enough to accomplish the finishing and polishing of the composite restoration, saving time and abrasive materials during this operation
- High resistance tips that can be used several times;
- Available in 3 shapes
- The mandrel and tips are autoclavable



Praxis

- For composite restoration finishing and polishing
- Pop on mandrel system, easy and fast for the installation of the discs as well as the inversion of the abrasive surface to mesial and distal positions
- Ultra-fine discs very flexible and versatile
- Four grits codified by colors, in sizes of 1/2' and 3/8'




CONSERVATIVE & ENDODONTICS



CONFIDENT
SALES INDIA PVT LTD


Superfix

- The Superfix system is used for resin restoration finishing and polishing
- Versaplus mandrel system that works by pressure, it does not have metallic core, what permits to use all the abrasive surface
- Ultra-fine discs, very flexible and versatile
- Four grits codified by colors, in sizes of 9,5mm and 12,7mm



TDV sandpaper disc

- These discs are indicated for dental finishing and polishing.
- Ultra-fine discs, very flexible and versatile
- Excellent cost-benefit relation
- Four grits codified by colors, in sizes of 12mm, 16mm and 19mm



Felt disc

- The felt discs are indicated for polishing with dental polishing paste.
- Perfect brightness on all type of surfaces
- Applicable with all sort of polishing pastes
- Available in sizes of 8mm and 12mm





CONFIDENT
SALES INDIA PVT LTD

Impregnated felt disc

- Natural felt discs with a 12 mm diameter, abrasive and lubricant impregnated, for compound, amalgam, and metal polishing
- Superb final polishing of composites and shiny surfaces
- Impregnated with ultra-fine abrasive dehydrated polishing paste; activated with a few drops of natural water
- Available in 12mm size



Polimax

- Superb final polishing of composites and shiny surfaces like acrylics, amalgam, metal and dental enamel.
- Impregnated with ultra-fine abrasive dehydrated polishing paste;
- Activated with a few drops of natural water.



Polifix

- Polifix discs are used for restoration material polishing with aid of a polishing paste.
- Perfect brightness on all type of surfaces
- Applicable with all type of polishing pastes



info@confidentonline.com

Toll Free : 1800-425-3366



Composite preparation:

Fotoplus:

- Developed to prevent a marginal leakage on resins Class II restorations providing good marginal sealing
- Translucent instruments that allows to light cure through their tips while the composites are being hold under pressure
- 2 instruments with 4 different tips


Esthetic plus:

- Anti-adherent silicone brushes for composites sculpting
- Hard and soft tips to find the best compatibility with the hardness of the composite .
- Tips available in four different shapes

Contact Gold:

Contact+ Gold is a unique instrument conceived to produce consistent contact points with ideal shape and position in class II restorations with composite resins.

- Translucent instruments that allows to light cure through their tips while the matrix is being hold attached to the adjacent tooth
- One instrument for large cavities (L) and one for small cavities (S)





CONFIDENT
SALES INDIA PVT LTD

Mandrel:

Metallic/Plastic Versaplug mandrel

Metallic mandrel to be used with the products Optimize, Superfix, Polifix, Polimax and Profilax.



Praxis mandrel

Mandril metal type **Soflex Pop-on***, in steel, to be used with Praxis discs.



Regular Mandrel:

Disk can be mounted on and demounted from the mandrel easy. These are of stainless steel with long 9mm screw and little washer.

Available in 2 types:

- Straight handpiece mandrel
- Contrangle handpiece mandrel.

