

 <p><b>A INSTRUMENT</b> For removing supragingival deposits in all quadrants</p> <p><b>DS-001A</b></p> <p><b>Rs : 3,400/- + Tax</b></p>	 <p><b>B INSTRUMENT</b> For cleaning lingual surfaces with heavy stain and calculus</p> <p><b>DS-003A</b></p> <p><b>Rs : 3,400/- + Tax</b></p>	 <p><b>C INSTRUMENT</b> For removing deposits and stain – especially suitable for removing bonding residue e.g. in orthodontics</p> <p><b>DS-002A</b></p> <p><b>Rs: 4,800/- + Tax</b></p>	 <p><b>P INSTRUMENT</b> For removing supra and subgingival calculus in all quadrants</p> <p><b>DS-011A</b></p> <p><b>Rs : 4,800/- + Tax</b></p>
<p><b>DS-016A</b></p>  <p><b>Rs: 4,800/- + Tax</b></p> <p><b>PS INSTRUMENT</b> For removing calculus in deep periodontal pockets in all quadrants “EMS PERIO SLIM TIP HAD BEST INTERPROXIMAL AND SUBGINGIVAL ACCESS” – CRA (Clinical Research Associates) Newsletter, USA, June 1998</p>	 <p><b>RT1 INSTRUMENT</b> Diamond-tipped instrument for locating canals and removing calcifications in the coronal third of the root canal system</p> <p><b>DS-066A</b></p> <p><b>Rs: 7,000/- + Tax</b></p>	 <p><b>RT2 INSTRUMENT</b> Diamond-tipped instrument for removing ledges. Creates straight canal for easy viewing of separated instrument</p> <p><b>DS-067A</b></p> <p><b>Rs : 7,000/- + Tax</b></p>	 <p><b>RT3 INSTRUMENT</b> Very long and fine instrument for the removal of broken instrument tips. Use with 180° file holder</p> <p><b>DS-062</b></p> <p><b>Rs : 7,500/- + Tax</b></p>
 <p><b>120° FILE HOLDER</b> For the flexibility to use different endo files, specially suited for anterior and premolar areas. Also suitable for RE2 and Berutti instruments</p> <p><b>DS-010</b></p> <p><b>Rs.: 7,000/- + Tax</b></p>	 <p><b>90° FILE HOLDER</b> For the flexibility to use different endo files, specially suited for the molar area</p> <p><b>DS-012</b></p> <p><b>Rs: 7,000/- + Tax</b></p>	 <p><b>180° FILE HOLDER</b> For use with RT3, Berutti and RE2 instruments</p> <p><b>DS-069</b></p> <p><b>Rs: 7,000/- + tax</b></p>	 <p><b>ESI INSTRUMENT</b> Non-cutting NiTi Non-cutting NiTi instrument for the effective cleaning and irrigation of the root canal system</p> <p><b>FS-279</b></p> <p><b>Rs: 7,500/- + tax</b></p>
 <p><b>D INSTRUMENT</b> For loosening crowns, bridges and pins placed previously with phosphate cement</p> <p><b>DS-005A</b></p> <p><b>Rs. 5,500/- + Tax</b></p>	 <p><b>H INSTRUMENT</b> For the lateral thermal condensation of gutta-percha</p> <p><b>DS-004A</b></p> <p><b>Rs. 5,500/- +Tax</b></p>	 <p><b>FLAT WRENCH</b> To change instruments for all three file holders</p> <p><b>DT-018</b></p> <p><b>Rs. 1,100/- + Tax</b></p>	 <p><b>RE2 INSTRUMENT</b> Instrument for retrograde root-canal isthmus preparation. The surface is coated with titanium nitride to improve abrasion. Use with 120° or 180° file holders</p> <p><b>DS-061</b></p> <p><b>Rs. 8,000/- + Tax</b></p>

 <p><b>BERUTTI INSTRUMENT</b> Diamond-tipped instrument for the retrograde preparation of root canals. A tapered design allows minimally invasive access and therefore fewer postoperative complications. Use with 120° or 180° file holders</p> <p><b>DT-060</b> Rs. 7,500/- + Tax</p>	 <p><b>STEEL FILES/K-FILES/ ENDO FILES</b> K-Files for cleaning and disinfecting the root canal system [ &gt; Use files in 90° or 120° file holders] 6 endo files ISO 15 (white) DT_006 6 endo files ISO 20 (yellow) DT_007 6 endo files ISO 25 (red) DT_008 6 endo files ISO 30 (blue) DT_009 6 endo files ISO 35 (green) DT_010 1 set endo files ISO 15/20/25/30/35 mm, 6 pcs each</p> <p><b>DT-011</b> Rs. 5,500/- + Tax</p>	 <p><b>PL1 INSTRUMENT</b> Left-curved instrument for periodontal debridement – especially suitable for difficult-to-reach interproximal surfaces</p> <p><b>DS-030A</b> Rs. 5,100/- + Tax</p>	 <p><b>PL2 INSTRUMENT</b> Right-curved. Otherwise same as PL1 instrument</p> <p><b>DS-031A</b> Rs. 5,100/- + Tax</p>
 <p><b>PL3 INSTRUMENT</b> Instrument for cleaning and disinfecting periodontal pockets – especially suitable for follow-up and recall visits</p> <p><b>DS-032A</b> Rs. 5,100/- + Tax</p>	 <p><b>PL4 INSTRUMENT</b> Left-curved instrument with small ball [ &gt; 0.6-mm diameter ] on instrument tip for the periodontal debridement of furcations and concavities</p> <p><b>DS-033A</b> Rs. 5,100/- + Tax</p>	 <p><b>PL5 INSTRUMENT</b> Right-curved. Otherwise PL4 instrument</p> <p><b>DS-034A</b> Rs. 5,100/- + Tax</p>	 <p><b>PF INSTRUMENT</b> Diamond-tipped instrument for enlarging f issues</p> <p><b>DS-062A</b> Rs. 6800/- + Tax</p>
 <p><b>FS-295</b> M.R.P : 10,500/-</p> <p><b>PI INSTRUMENT</b> For the effective and gentle cleaning of implants and corresponding restorations</p> <p><b>STARTER KIT</b> 4 PI instruments 1 f file holder (120°) 1 f lat wrench</p> <p><b>REFILL</b> 4 PI instruments</p> <p><b>DS-065A</b> Rs : 7,200/- +Tax</p>	 <p><b>HPL3/DPL3 INSTRUMENTS</b> Diamond-tipped perio instruments for in-depth root debridement using direct vision [ &gt; "Flap" ] – for smoothing restoration overhangs as well as expanding root furcations</p> <p><b>DS-065A</b> Rs. 9,000/- + Tax</p>	 <p><b>SP INSTRUMENT</b> For bonding ceramic inlays, onlays and veneers [ &gt; without SP cap ]</p> <p><b>DS-015A</b> Rs. 5,800/- + Tax</p>	